



## **How to use this CD**

### **Step 1**

#### **Install the USB Drivers**

Connect your monitor or wireless receiver to the computer using the supplied USB cable. The PC will detect the new hardware and what happens next depends upon the operating system you are running under. Follow the document "Lightning First Setup" on the CD or download it from the website at <http://www.vetronic.co.uk/products/article/lightning>.

Any drivers that are required can be found on your CD in the FTDI Drivers section under the appropriate operating system folder.

### **Step 2**

Follow the section in the Lightning First Setup CD to install the Vital Monitor and VitalStore software programs