



ICON



In Circuit Oxygen aNalyser

Want to know the oxygen concentration
in your Oxygen Tent?

ICON is your solution



- Simple to use
- Connect with 22mm and 15mm standard connectors
- Use with Oxygen concentrators
- Measure Fresh Gas Outlet concentrations
- Monitor Oxygen Tent concentrations
- Very long battery life
- Very long sensor life
- Simple to re-calibrate

Icon is a new in-circuit oxygen monitor designed for use in veterinary operating theatres and intensive care wards.

The Icon unit can be placed on the Fresh Gas Outlet of an anaesthetic machine or within the rebreathing limbs of a circle system, wherever the oxygen concentration needs to be known. Icon does not need a flow of gas to function and can be used inside an oxygen tent to monitor oxygen levels for high-dependency patients.

The Icon unit uses a single AA battery that has a lifetime of approximately 1 year when used for 8 hours a day every day. The sensing element is a long-life galvanic sensor that has a working life of approximately 2 years.

£275.00