

Impact-III FAQ

Q: What sampling method is used by the Impact-III capnography module?

A: The Impact-III capnography module uses a sidestream method for CO2 analysis

Q: What is the sidestream sampling rate of the Impact-III capnography module?

A: Fixed at 50ml/minute for all species

Q: What animal range is supported by the Capnography module?

A: There is no upper weight limit. The lower weight limit is about 150g but this does rely on a low dead-space ET Connector and careful attention to additional dead space.

Q: Do I need to buy all 3 modules that Impact-III supports?

A: No. It is possible to buy Impact-III in the following combinations: Pulse-Ox+CO2+Temperature, Pulse-Ox + Temperature, CO2 only, CO2 + Temperature

Q: Is the Impact-III battery powered?

A: Yes. There are internal NimH batteries so the unit can be run either from the batteries or from an attached mains charger

Q: How long will the unit operate when running from batteries?

A: With all modules fitted, the unit will run for up to 6 hours. This will increase by an hour or two if either the CO2 or Pulse-Ox modules are not fitted.

Q: Does the Impact-III have a touch screen?

A: No. All features and menu options are accessed using a simple 2-button keypad.

Q: Can the Impact-III be connected to a PC?

A: Yes, using an optional wireless link. A direct USB connection is not available as the unit has been designed for mainly portable use

Q: Where can I get the accessories I need for Impact-III?

A: We stock and supply all of the accessories for Impact-III in the UK and supply all of our distributors, so these will be available either directly from the UK or through the supplier of your Impact-III

Q: Do you repair the item in the UK?

A: Yes. All repairs are carried out in the UK with the exception of custom modules such as the Phasein Gas Analysers which we send to Sweden for repair

Q: Is it designed specifically for animals?

A: Yes. All of our equipment is designed specifically for the animal market. We do not adapt medical devices for veterinary use.

Q: What warranty do you offer?

A: We offer a standard 1 year warranty covering parts and labour. Ask us for more details.